

PVI00587EN_V01

Single-Chamber Sampler EK 100

Space-saving design consisting of drive unit and release chamber

- Compact and space-saving single-chamber design
- Suitable for different product flow properties
- Manual or automatic sampling
- ~15-25 cm³ per sample (15-25ml)
- Screw thread for custom sample bottles
- Endless-Foil unit for containment applications
- Easy disassembly and easy cleaning



Function endless foil unit

The system operates in closed mode. After having reached the total sample quantity, the filled bag is pulled downwards a little and closed with a clip. In order to separate the bag and have the next bag ready for the next sampling, another clip is placed underneath the first and then it is possible to cut between the two clips and take tube the filled bag.

